BORT OmoControl

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Med
Mitten im
gesunden
Leben

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REF 120 700, 120 500/600

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	03
EN	English	Instructions for use	12
FR	Français	Mode d'emploi	21
ES	Español	Instrucciones de uso	30
IT	Italiano	Istruzioni per l'uso	39
NL	Nederlands	Gebruiksaanwijzing	48

PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a shoulder joint brace for the physiological guidance and/or stabilisation of the shoulder joint. Joint-damaging or dislocation-promoting movements of the upper arm in the shoulder joint should be hindered to a great extent.

Indications

Flaccid and painful movement restrictions of the shoulder and arm muscles after: slight to moderate instability of the shoulder joint, conservative in case of permanent dislocation of the shoulder, inflammation and degeneration

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- do not close it too tightly during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions: Information for the patient

The instructions describe putting it on when sitting

1
Open the BOA® pulley system by pulling the rotary knob (clicking noise). Pull the two brace elements apart until the entire rope is pulled out of the BOA® system. For simplified handling, one of the brace elements can be fixed beneath your thigh at the same time. Then, press the BOA® system's rotational mechanism back together again – when engaging, you can perceive a clicking sound – the rope take-up is now secure during the later course.

Attachment process

2
Tip: The arm cuff can be fixed beneath the thigh. Pull the Velcro fasteners fully out of the arm cuff guide loops.

3
Tip: Fix the long end of the chest strap beneath the thigh. Guide the index finger through the finger flap and pull the guide loop up to the twist and click fastener using the thumb and middle finger. This extends the chest strap to its full length. The twist and click fastener remains closed.

4
Loosen the axillary strap, this enables easier passage of the arm affected during the later course.

5
Open the blue quick-release fastener by releasing it. The rear grey quick-release fastener remains closed.

6
Guide the arm affected below the chest strap and through the axillary strap.

7
Use the healthy side of the arm to grasp the upper edge of the shoulder cuff and pull the shoulder cap upwards over the shoulder. Here, the chest strap is also guided upwards.

8
Guide the chest strap over your head and over your healthy arm. When pulling it over, pay attention not to twist the shoulder strap.

9
Tighten the chest strap pleasantly tight and attach it using the Velcro fastener.

10 Ensure here that the upper guide of the pulley system is positioned on the roof of the shoulder (acromion). This is an important basic position for guaranteeing optimum functionality.

11 Close the axillary strap. Ensure a firm, but not too tight fit.

12 Clamp the blue arm pad quick-release fastener into the lower slit of the guide loop.

13 Then, position the elbow in the arm pad. Here, the arm pad can be placed on the thigh. For self-control, the elbow should be visible in the windowed area of the pad. After you have guided the Velcro fasteners through the guide loops, they can be closed with a firm, pleasant pull.

14

Release the click fastener previously fixed to the axillary strap guide loop and join it to the counterpart on the arm pad fastener by clasping it. The shoulder pad pulley system is then connected to the arm pad.

Information: The matching click fasteners have the same colours. The outer fastener on the reverse side is grey, the inner close to the body blue.

15

Assume an upright position carefully. By turning the BOA® system clockwise, the arm can now be pulled towards the acromioclavicular joint.

16a

The fitting process is completed.

16b

Tip: You can find our video, which shows the attachment process, both on our homepage and via the product-specific QR code.

Removal

To remove the OmoControl, first open the twist and click fastener, then the attached blue click fastener, the Velcro fasteners and the axillary strap. Then, slip the OmoControl off on the side affected with the healthy arm.

Fitting instructions: Information for the technician

1

Open the BOA® pulley system by pulling the rotary knob (clicking noise). Pull the two brace elements apart until the entire rope is pulled out of the BOA® system.

2

Remove the quick-release fasteners which connect the shoulder pad and the arm pad.

3

Open the twist and click fastener on the chest strap.

4

Pull the axillary strap to its full length and attach it again firmly.

5

Guide the affected side of the arm through the axillary strap and pull the shoulder pad onto the roof of the shoulder. Pay attention to flush fitting placement on the shoulder.

6

Close the chest strap beneath the axilla via the twist and click fastener and pull the belt tight using the thumb loop. The thumb loop can be attached to the chest strap using the Velcro surface on the thumb loop.

7

If required, the length of the chest strap can be shortened via the crocodile fastener, thus avoiding superfluous strap length in the thumb loop area.

8

Close the axillary strap. It can be shortened via the crocodile fastener. When closing it, pay attention to a firm, but not too tight fit.

9a

Then, position the elbow in the arm pad. Adapt the fasteners on the arm pad using the crocodile fastener. For self-control, the elbow should be visible in the windowed area of the pad. The Velcro fasteners can then be closed with a firm, pleasant pull.

9b

10a

10b

Join the colour-marked click fasteners with their counterpart fasteners on the arm pad. The outer fastener on the reverse side is grey, the inner close to the body blue.

11

Finally, engage the BOA® wheel again by pressing it. By turning the BOA® system clockwise, the arm can now be pulled towards the acromioclavicular joint.

12a

12b

The fitting process is completed.

Note: The arm pad and shoulder pad can be ordered separately as modules and have the following article numbers: arm pad [REF](#) 120 500 (module A) and shoulder pad [REF](#) 120 600 (module B). You can find the appropriate size selection in the size chart in the programme booklet or the price list.

Material composition

51% polyamide (PA), 32% polyester (PES), 10% viscose (CV), 3% elastane (EL), 3% polyurethane (PUR), 1% cotton (CO)

Cleaning information

 Delicate wash 30 °C  Do not bleach (only wash with colour or mild detergents)

 Do not dry in a tumble dryer  Do not iron  Do not clean chemically

Close the Velcro fastener to avoid damage to other laundry items. Do not use fabric conditioner.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: <http://www.bort.com/conformity>

Stand: 11.2023

 Medical device |  Single patient – multiple use

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