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
## BORT Asymmetric Sport with a grip top

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### **Intended purpose**

This medical device is a textile brace for facilitating movement of the patella. The pressure pad on the kneecap and the adjustable tension element can achieve correction and relief of the kneecap.

### **Indications**

Patellofemoral pain syndrome, high riding patella, patella lateralisation, condition after lateral release OP, recurrent or permanent patella instability, meniscopathy.

### **Contraindications**

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

### **Application risks/Important notes** ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds

- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

### Fitting instructions

- 1** Fold the grip top outwards, grasp the upper edge of the support.
- 2** It's best to put the aid on when sitting. Slip the support over the forward section of the foot and then pull it up to the knee.
- 3** Position the support on the knee joint so that the kneecap (patella) is supported at the side by the integrated silicone half-ring.
- 4** In the next step, remove the upper Velcro surface restraint and re-attach it with Velcro in the direction of the inside leg under slight tension. Then attach the lower restraint in the direction of the leg under slight tension.
- 5** Then, flap the grip top back towards the leg so that the silicone nodules are directly on the skin.

### Removal

- 6** To remove the support, please fold the grip top outwards
- 7** Then, remove the correction constraint, grasp the lower edge of the support at the side on the splints and pull it downwards over the foot. Close the Velcro fastener in order to avoid damaging other items, e.g. to further clothing.

## Material composition

Polyamide (PA), elastane (EL), silicone (SI)

The sewn in textile label provides the precise material composition.

## Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Close the Velcro fastenings before washing to avoid damage to other laundry items.

## Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

**Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

**Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link:  
[www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2024

**MD** Medical device |  Single patient multiple use

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