
BORT MyoActive Sport Wadenstütze

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
Sport

Illustrationen
Figures


1


2



3

4

5

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	04
EN	English	Instructions for use	08
FR	Français	Mode d'emploi	12
ES	Español	Instrucciones de uso	16
IT	Italiano	Istruzioni per l'uso	20
NL	Nederlands	Gebruiksaanwijzing	24
CS	Čeština	Návod k použití	28
ET	Eesti	Kasutusjuhend	32
PL	Polski	Instrukcja użytkowania	36
RO	Română	Instrucțiunile de utilizare	40

PDF: ga.bort.com

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BORT MyoActive Sport Calf Support

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a lower leg support for the compression of soft tissue with diagonal and elastic restraints as well as distal and proximal anchors.

Indications

For the compression of soft tissue on the lower leg in case of acute ruptures of the muscle fibre and/or acute muscular injuries

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used

- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

1 Open both diagonal cross restraints and let the restraints hang downwards. Slip the support with the larger of the two openings (the top) over the foot first. Put the support on when sitting. Here, the relief zone / function zone (F) should be at the height of the injured region

Information: The diagonal cross restraint reinforcement is in the calf region. The front support seam should run centrally along the edge of the calf.

2 On the inner side of the cross restraint, you can see the green inner velour surface with the directional arrows and marking strips as a so-called “pull stop”. Then, pull the first cross restraint upwards to the desired marking strip (in the direction of the arrow). Do not exceed the maximum tension shown by the widest, red marking strip here.

3 Then, place the first cross restraint (in the guide loop) and attach it diagonally to the lower velour surface using Velcro.

Information: If the velour restraint is too long, you can shorten it by removing the crocodile fastener. Then, attach the crocodile fastener at the end of the velour.

4 Repeat the process with the second cross restraint. Here, the restraints cross as an X. It makes no difference which cross restraint you attach first.

5 If required, the circumferential and elastic belt fasteners provided can be attached to the velour anchors under tensile stress. Ensure that the longer restraints are attached above and the shorter restraints below.

Removal

To remove, if used, loosen the elastic circumferential belt fasteners (anchors) if used. Open the cross fasteners and take hold of the knitted fabric on the support to pull it downwards over your foot.

Material composition

Cotton (CO), polyamide (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

Cleaning information

Delicate wash 30 °C Do not bleach (only wash with colour or mild detergents) Do not dry in a tumble dryer Do not iron
 Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fasteners to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

MD Medical device |  Single patient – multiple use

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