



# BORT SELLADUR<sup>®</sup>

## THUMB ORTHOSES



Injection molded  
models in a unique  
perforated design

Made in Germany

# COMFORTABLE, UNIVERSAL

## THANKS TO PERFORATION AND SPECIAL SHAPING

Pains have many causes and are often unavoidable. Immobilization supports the healing process, and it frequently also brings relief of pain. For this purpose we developed orthoses.

BORT SellaDur® injection molded thumb orthoses are now more comfortable to wear. And these are the decisive advantages:

- > Very comfortable to wear due to perforation, which allows plenty of air to reach the covered areas of skin. This provides optimal ventilation.
- > Due to the special molding material the SellaDur® cannot break when being fitted or if inadvertently squeezed.
- > In the critical thumb-forefinger-area special high formed edge avoids pressure spots.
- > Washable and hygienic.

### Indications:


Immobilization of carpal and metacarpal joint, ulnar ligament lesions (ski thumb), partial immobilization, rhizo-arthritis.

### Sizes:

BORT SellaDur® thumb orthoses are available in three different sizes.

### Measuring point:

Width of hand at metacarpal knuckle joint.



Sizes	Hand Width in cm
small .....	7,0 – 8,0 cm
medium.....	8,0 – 9,0 cm
large .....	9,0 – 10,0 cm

Please indicate left or right.



Art. no. 105 400

**BORT Sella Dur®**

Standard version. Color: blue.



Art. no. 105 600

**BORT Sella Dur® FLOC**

Additional flocking with textiles provides optimal wearing comfort. Color: grey.



**BORT GmbH**

P.O. Box 1330  
D-71367 Weinstadt  
medical@bort.de  
www.bort.de

QUALITY  
MANAGEMENT  
certified acc. to  
medical product guidelines  
93/42/EWG



QUALITY  
MANAGEMENT  
we are certified  
regular voluntary  
control acc. to ISO 9001:2000

