

NOVELTIES 2025



Pain relief in case of overstress

BORT SellaTex® Plus - Quervain

Wrist brace for immobilisation and stabilisation in case of Tendovaginitis de Quervain



With medial opening

Med

BORT DigiSoft® Vario Finger Brace

Sport

BORT MyoActive Sport Calf Support and Asymmetric® Sport Knee Support with a Grip Top

StabiloGen® Eco

Our bestseller with a Grip Top Short
Now available in black

BORT. Benefit at your side.



NEW DESIGN

BORT SellaTex® Plus - Quervain

Wrist brace with medial opening for immobilisation and stabilisation



REF 112 160

PRODUCT-BENEFITS

- Wrist brace for stabilising the metacarpophalangeal joint, saddle joint and end joint of the thumb
- Palm and dorsal splints for wrist joint and metacarpus stabilisation
- **Can be opened fully**
- Length: approx. 20 cm

COLOUR: ■ grey

THE MED-PLUS

- **Large opening in the region of the problematic medial area** in case of Tendovaginitis de Quervain
- Easy downgrade option for the thumb splint within the course of therapy
- Long thumb element can be individually positioned, can be shaped, fully removable
- Pleasant wearing comfort



Thumb splint which can be downgraded for easy adaption during the course of therapy

INDICATIONS

Tendovaginitis de Quervain, pre-operative, post-traumatic, post-operative, e.g. after thumb ligament reconstruction, immobilisation for rheumatic inflammatory processes, rhizarthrosis and for degenerative diseases, CTS

SIZES

Circumference of wrist


	x-small	_____	up to 15 cm
	small	_____	15 – 17 cm
	medium	_____	17 – 19 cm
	large	_____	19 – 21 cm
	x-large	_____	21 – 23 cm

please indicate left or right

Special wrist brace for conservative Tendovaginitis de Quervain therapy

OVERSTRESS & PAIN

Carpal tunnel syndrome in the region of the thumb side of the wrist is a chronic overstress symptom and is primarily caused by frequent abduction (movement of the thumb away from the palm).



IMMOBILISATION & RELIEF

For conservative therapy, the hindrance of overburdening and immobilisation of the joint region are focussed. The **BORT SellaTex® Plus – Quervain** Wrist Brace has a stabilising effect and immobilises the thumb joint. As a result, it can hinder movements which trigger pain.



NEW

BORT DigiSoft® Vario Finger Brace

For reliable immobilisation



REF 112 790


PRODUCT-BENEFITS

- Finger brace for immobilisation of the interphalangeal joints
- Aluminium splint, adjustable with soft padding
- Velcro fasteners, variably positionable
- **Can be opened fully**

COLOUR: ■ grey

THE MED-PLUS

- Modular brace system:



Wrist stabilisation, downgradeable

SIZES Metacarpus circumference

	Size 1	17 – 19 cm
	Size 2	19 – 22 cm
	Size 3	22 – 26 cm

right and left symmetrical

Med

NEW DESIGN

BORT StabiloGen® Eco with a Grip Top Short


Now finally available in black



REF 114 520 HK

PRODUCT-BENEFITS

- Classical knee support for compression of soft tissue with a visco-elastic pad for the stabilisation and relief of the knee joint
- **With additional silicone grip top in all sizes**
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- One lateral spiral spring stay on each side
- **Double Plus size system, therefore an excellent fit**

COLOURS: ■ silver ■ black 


THE ECO-PLUS

- **Breathable functional comfort zone in the popliteal fossa**
- Silky, soft surface structure
- **Short version** in the thigh area, no extension thanks to the silicone grip top
- Pleasant wearing comfort
- Length of the support in size 3: 32 cm

SIZES d-measure f-measure

	Size 1	26 – 29 cm	32 – 36 cm
	Size 2	29 – 33 cm	36 – 40 cm
	Size 2+	29 – 33 cm	40 – 44 cm
	Size 3	33 – 37 cm	40 – 44 cm
	Size 3+	33 – 37 cm	44 – 48 cm
	Size 4	37 – 41 cm	44 – 48 cm
	Size 4+	37 – 41 cm	48 – 52 cm
	Size 5	41 – 45 cm	52 – 56 cm
	Size 5+	41 – 45 cm	56 – 60 cm
	Size 5++	41 – 45 cm	60 – 64 cm
	Size 6	45 – 49 cm	60 – 64 cm
	Size 6+	45 – 49 cm	64 – 68 cm
	Size 6++	45 – 49 cm	68 – 72 cm
	Size 7	49 – 53 cm	68 – 72 cm
	Size 7+	49 – 53 cm	72 – 76 cm

right and left symmetrical



Med

NEW DESIGN

Product modifications

Colour changes



BORT MalleoStabil® Brace / MalleoStabil® SOFT Brace / MalleoStabil® Air/Gel Brace

REF 100 500 / 100 520 / 100 540

- New colour: ■ grey

Design changes



BORT StabiloPro® Knee Brace

REF 182 150

- With a new polar curve joint based on the knee joint's physiological rolling and sliding movement



BORT StabiloPro® Knee Brace Open Style

REF 182 300

- With a new polar curve joint based on the knee joint's physiological rolling and sliding movement



BORT Generation Knee Brace

REF 215 400

- With a new polar curve joint based on the knee joint's physiological rolling and sliding movement



BORT Osgood-Schlatter Knee Support

REF 114 510

- Now also available in sizes 6 and 6+



BORT CoxaPro Hip Brace

REF 201 400

- New with various f-measures



BORT TalofX®

REF 054 800

- New style
- New Velcro fastener system
- Can be opened fully for easy fitting

NEW

BORT MyoActive Sport Calf Support

Active support for ruptures of the muscle fibre and muscular injuries



REF 116 500 SP

PRODUCT-BENEFITS

- Lower leg active support with circular straps and cross restraints
- Anatomically contoured knit
- Relief zones in the injured region
- Lengths: 24 cm (small), 25 cm (medium), 26 cm (large; x-large)

COLOUR: ■ black/green

THE SPORT-PLUS

- Mechanism of action from compression, conventional tape and kinesio tape
- Proprioceptive effect
- Modern sporty design
- Latex-free fabric

Sport

SIZES		Circumference of calf
	small	up to 33 cm
	medium	33 – 38 cm
	large	38 – 43 cm
	x-large	43 – 47 cm

right and left symmetrical

NEW

BORT Asymmetric® Sport Knee Support with a Grip Top

For correction of the patella glide path, also available in Plus and Double Plus sizes



REF 114 700 SP H

PRODUCT-BENEFITS

- Knee support for influencing the patella glide path
- SOFTflex 3D knitted fabric for more elasticity, stability and very comfortable to wear
- With additional silicone grip top in all sizes
- Double Plus size system, therefore an excellent fit
- One lateral spiral spring stay on each side
- Compression class II

COLOUR: ■ black/green

THE SPORT-PLUS

- Pad pressure individually and directly adjustable using two Velcro straps
- COOLMAX®: Wicks moisture away from the body
- Ideal fit for sports activities
- Latex-free fabric
- Modern sporty design (with a reflecting label)
- Support length in size 3: 35 cm

SIZES		d-measure	f-measure
	Size 1	26 – 29 cm	32 – 36 cm
	Size 2	29 – 33 cm	36 – 40 cm
	Size 2+	29 – 33 cm	40 – 44 cm
	Size 3	33 – 37 cm	40 – 44 cm
	Size 3+	33 – 37 cm	44 – 48 cm
	Size 4	37 – 41 cm	44 – 48 cm
	Size 4+	37 – 41 cm	48 – 52 cm
	Size 5	41 – 45 cm	52 – 56 cm
	Size 5+	41 – 45 cm	56 – 60 cm
	Size 5++	41 – 45 cm	60 – 64 cm
	Size 6	45 – 49 cm	60 – 64 cm
	Size 6+	45 – 49 cm	64 – 68 cm
	Size 6++	45 – 49 cm	68 – 72 cm
	Size 7	49 – 53 cm	68 – 72 cm
Size 7+	49 – 53 cm	72 – 76 cm	

please indicate left or right



COOLMAX® is a trademark of The Lycra Company.

Sport

NEW DESIGN

BORT X-Walker

The brace programme for the care of foot and ankle joint injuries



BORT X-Walker long

REF 101 300




BORT X-Walker short

REF 101 350



BORT X-Walker Achillo

REF 101 320



Extensive vacuum system

for flush fitting integration and stabilisation. The vacuum cushion integrated in the padding liner filled with granulate can be perfectly adapted to the shape of the foot or lower leg using the vacuum pump provided. Pressure sores can be avoided due to the even and individually adjustable bedding.



Spare parts for the BORT X-Walker



Liner-Set
(liner cover, vacuum cushion and pump)



Liner cover



Vacuum pump



Heel wedges set



Vacuum cushion

⊕ Accessories for all BORT X-Walkers

REF 100 370
BORT Understocking



REF 100 380
BORT Protective cover against wetness



REF 100 390
BORT Protective cover for night-time



Med

All information without guarantee, error and changes reserved. Status: January 2025.

BORT. Benefit at your side.

BORT GmbH | P.O. Box 1330 | D-71367 Weinstadt | Am Schweizerbach 1 | D-71384 Weinstadt
Phone +49 (0) 7151 99200-0 | Fax +49 (0) 7151 99200-50 | www.bort.com | medical@bort.com