

The new BORT StabiloGen Eco

For your standard treatment in the knee area



Med

In the midst
of healthy
life



BORT. Benefit at your side.



INDIVIDUAL FIT

The even degree of extensibility of the knitted compression fabric adapts to the body shape.

FUNCTIONAL EFFECT

The compressive, massaging effect of the knitted fabric provides permanent proprioceptive, neuromuscular stimuli. This means that muscles are activated, the knee joint can be stimulated and pain relief be provided (gate control).

VISCOELELASTIC PAD

For stabilisation and relief of the knee joint

HIGH DIMENSIONAL STABILITY

With a lateral spiral stay on each side for constant fit of the support on the leg.

BORT StabiloGen **Eco** - the support for your standard treatment

Versatile functions and product variants

PRODUCT-BENEFIT

- Classical knee support for compression of soft tissue with a viscoelastic pad for stabilisation and relief of the knee joint
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- One lateral spiral stay on each side
- Plus-size system results in a good fit
- Compression class II

COLOURS: ■ skin-tone ■ silver ■ black

THE ECO-BENEFIT

- Breathable functional comfort zone in the popliteal fossa
- Velvety soft surface structure
- Silicone adhesive border from size 5

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissue, recurrent articular effusion, e.g. in case of gonarthrosis, arthritis, ligament instability, meniscopathy.

COMFORT ZONE POPLITEAL FOSSA

Breathable functional comfort zone.



SOFTflex knitted fabric: MODERN KNITTING TECHNOLOGY WITH A FEEL-GOOD EFFECT

Velvety soft surface structure – which is comfortable to wear on the skin.



The original: Plus-size system from BORT

Better fit, better function, more comfortable to wear

BORT PLUS-SIZE SYSTEM FOR THE THIGH AREA

- 5 Plus-sizes – As intermediate sizes in f-measure from size 2+ (40 cm) to size 6+ (68 cm) for a considerably better fit in the thigh area.
- Better function and more comfortable to wear.

SILICONE ADHESIVE BORDER FROM SIZE 5



The multi-variant solution in the series portfolio

BORT StabiloGen Eco

REF 114 520



PRODUCT-BENEFIT

- The classic product in standard format for the large majority of (straight) leg shapes
- Length of the support in size 3: 32 cm

BORT StabiloGen Eco
with a Grip Top Short

REF 114 520 HK

NEW



Additional silicone grip top in all sizes,
height 5 cm

PRODUCT-BENEFIT

- Extensive thigh frame for short thigh lengths
- Additional silicone grip top, height 5 cm, in all sizes
- Length of the support in size 3: 32 cm

BORT StabiloGen Eco

with a Grip Top

REF 114 520 H



PRODUCT-BENEFIT

- Extensive thigh frame
- Additional silicone grip top, height 5 cm, in all sizes
- Length of the support in size 3: 37 cm

Product variants and size overview



BORT StabuloGen Eco

REF 114 520

COLOURS:

■ skin-tone ■ silver ■ black



BORT StabuloGen Eco

with a Grip Top
REF 114 520 H

COLOUR: ■ silver



BORT StabuloGen Eco

with a Grip Top Short
REF 114 520 HK

COLOUR: ■ silver



BORT StabuloGen Eco

REF 114 520

COLOUR: ■ denim



BORT StabuloGen Eco

REF 114 520 KIDS

COLOUR: ■ blue

d-measure	f-measure	SIZES
23 – 26 cm	28 – 32 cm	Kids*(blue)
26 – 29 cm	32 – 36 cm	Gr. 1
29 – 33 cm	36 – 40 cm	Gr. 2
29 – 33 cm	40 – 44 cm	Gr. 2+
33 – 37 cm	40 – 44 cm	Gr. 3
33 – 37 cm	44 – 48 cm	Gr. 3+
37 – 41 cm	44 – 48 cm	Gr. 4
37 – 41 cm	48 – 52 cm	Gr. 4+
41 – 45 cm	52 – 56 cm	Gr. 5
41 – 45 cm	56 – 60 cm	Gr. 5+
45 – 49 cm	60 – 64 cm	Gr. 6
45 – 49 cm	64 – 68 cm	Gr. 6+

*Children's size only available for REF 114 520
right and left symmetrical
1 cm = 0,3937 inch



114 520

Measuring point:
d = 15 cm below the patella
f = 15 cm above the patella

BORT training app

10 selected exercises for beginners and advanced users to improve mobility and muscular stability



Off you go under:
uebungen.bort.com



MZ089/2020-11/001-GB

